



## **Autoclave TRBA-601**

[www.labdel.com](http://www.labdel.com) | [info@labdel.com](mailto:info@labdel.com)

## Autoclave TRBA-601

Benchtop Autoclave TRBA-601 is a pressure steam sterilizer, covered with a 304 hardened steel bowl and a manual door lock. It has a chamber volume of 16 L and an extreme working temperature of 136 °C. Having a security alert that demonstrates the fluctuations. In addition, the division, it is having presented three sterilizing plates. Handled with a fast cleansing of 15 to 30 min.

### Features:

Based on the sources, it is introduced with a wrapped and opened up system

Both systems having with the conclusion of 121 or 134°C

Sterilization is computer controlled

The chamber is made of extraordinary hardened steel

Working pressure time- 0.9 bar to 2.3 bar; greatest working temperature: 136°C

The division is presenting with three sterilizing plates

Security valve, manual door lock, fault different test and overheat protector

A fast sterilizing measure in 15-30 minutes

## Specifications :

<b>Appearance Size</b>	445×550×395 mm
<b>Frequency</b>	50/60 Hz ±1 Hz
<b>G·W</b>	32 kg
<b>Inside Dimension</b>	? 249×360 mm
<b>Package Dimension</b>	670×550×500 mm
<b>Power</b>	1800 W
<b>Tank Capacity</b>	16 Liters
<b>Voltage</b>	AC 220 V/110 V±10 %



Email: [info@labdel.com](mailto:info@labdel.com) | Website: [www.labdel.com](http://www.labdel.com)