



Centrifuge TRC-602

www.labdel.com | info@labdel.com

Centrifuge TRC-602

Centrifuge TRC-602 is a small device that can accommodate a variety of sample tubes. It has a maximum speed of 4000 rpm and an RCF of $1790\times g$. It has a rotor capacity of $20\text{ ml}\times 6$. It is very easy to use, comes with a switch to continue operating and set parameters, and is well suited for qualitative analysis.

Features:

It is having with a Scope of supply that includes the rotor

It can hold various sample tubes

It is having with a Small footprint

Very simple to operate

Having a Step speed adjustment

It Suits for qualitative analysis

Specifications :

Dimension	Ø250×220 mm
Maximum RCF	1790×g
Maximum RPM	4000 RPM, step speed adjustment
Rotor Capacity	20 ml×6
Timer	0 to 30min
Voltage	110V 60Hz /220V 50Hz
Weight	3.5 Kg



Email: info@labdel.com | Website: www.labdel.com