



## Centrifuge TRC-603

[www.labdel.com](http://www.labdel.com) | [info@labdel.com](mailto:info@labdel.com)

## Centrifuge TRC-603

Centrifuge TRC-603 with a powerful plastic sleeve and a steel lid, as well as a high brushless motor. It has a large LCD, a mechanical door lock, and a safety switch. It has a maximum rpm of 4000 and an RCF of 2325×g. It has an adjustable speed up and down and a standard rotor capacity of 20 ml ×12. It is appropriate for the qualitative analysis.

### Features:

It is a durable plastic container with a metal lid

Having a brushless high-speed motor

With an LED display, speed, and timer

It has over speed and imbalance warnings

It is built with a mechanical door lock and a safety switch are included

It displays the error code and allows you to adapt the speed up and down

It is appropriate for qualitative analysis

## Specifications :

<b>Dimension</b>	280×310×265 mm
<b>Maximum RCF</b>	2325×g
<b>Maximum RPM</b>	4000 RPM, step-less speed adjustment, LED display speed, and timer
<b>Rotor Capacity</b>	20ml ×12 (standard)
<b>Timer</b>	0 to 99 min
<b>Voltage</b>	110 V 60 Hz /220 V 50 Hz
<b>Weight</b>	9 kg



Email: [info@labdel.com](mailto:info@labdel.com) | Website: [www.labdel.com](http://www.labdel.com)