



Drying Oven TRDO-707

www.labdel.com | info@labdel.com

Drying Oven TRDO-707

Drying Oven TRDO-707 is a Microprocessor mechanism that has a self-tuning PID controller and a temperature range of 50 to 300 °C. It has a double glass window and a stainless-steel interior. The corners of the chamber are semicircles. A digital display setting and actual temperature, as well as a running timer ranging from 0 to 999 minutes, are included.

Features:

- It is with microprocessor control and a double-glazed window
- The real temperature and the temperature set just on the digital screen
- A timer is running on a digital display
- It is having an over-temperature alarm PID self-tuning
- Following the operation, an automatic alarm sounds
- The interior is stainless steel, semi-circular chamber corners

Specifications :

Dimension	720×780×750 mm
Power	1500 W
Temperature Range	50 to 300 °
Temperature Resolution	0.1 °
Temperature Fluctuation	±1 °
Timer	0 to 999 min
Voltage	110V 60Hz /220V 50 Hz
Weight	64 kg
Working Size	550×450×550 mm



Email: info@labdel.com | Website: www.labdel.com