



Digital Heating Water Bath TRHB-603

www.labdel.com | info@labdel.com

Digital Heating Water Bath TRHB-603

Digital Heating Water Bath TRHB-603 is a microprocessor controller with an inner digital display showing a temperature range of (RT+5) -300 °C. It holds with a capacity of 10 L. It is designed with a stainless steel interior with an aluminum ring cover.

Features:

It is equipped with a microprocessor control
mechanism with a stainless steel interior
Designed with a digital display inner temperature
It is covered with an aluminum ring
It is a thermostatic device

Specifications :

Capacity	10 L
Dimension	37x37x24 cm
Power	1000 W
Temperature Fluctuation	±1 °
Temperature Homogeneity	±1 °
Temperature Range	(RT+5) -300 °
Voltage	110V 60Hz /220V 50 Hz
Weight	5.5 kg
Working Size	30x30x12 cm



Email: info@labdel.com | Website: www.labdel.com