



Mini Horizontal Electrophoresis Cell TRMH-601

www.labdel.com | info@labdel.com

Mini Horizontal Electrophoresis Cell TRMH-601

Mini Horizontal Electrophoresis Cell TRMH-601 has primarily been used for rapid electrophoresis of small amounts of agarose gel containing DNA and RNA. It has the gel tray, lower body, upper body, gel-making kit, comb, and so on as the main components of the instrument. Built-in movable e-engine with no washer, simple maintenance, and no leakage. With the ability to separate as well as analyses up to 32 samples in 20 minutes, as well as a gel portion of 70×70 mm and 70×100 mm, it is filtered with high thermal stability and is non-deformed at 100°C.

Features:

In the composting of DNA molecules, it is beneficial to separate and analyze at most 32 samples in 20 minutes

Built-in movable e-engine with no washer, simple maintenance, and no leakage

At 100°C, the gel tray with a fluorescent measurement is resistant to heat and non-deformed

Gel irrigation can be performed before cooling the agarose

Specifications :

Buffering Liquid Volume	260 ml
Comb Specification	0.75mm: 9wells/16wells
	1.0mm: 9 wells/16 wells,
	1.5mm: 9 wells/16 wells
Gel area	70×70mm & 70×100mm
Weight	0.6 Kg



Email: info@labdel.com | Website: www.labdel.com