



NC Linear Shaker TRLS-601

www.labdel.com | info@labdel.com

NC Linear Shaker TRLS-601

NC Linear Shaker TRLS-601 has dual LCD screens for parameter setting and viewing. It can withstand a maximum load weight and operates in a timed and constantly functioning mode. The dual LCD screens separately display time and speed, with a speed range of 40 to 200 rpm and a time range of 1 minute to 19 hours 59 minutes.

Features:

There is a maximum load-bearing weight of 7.5kg allowed

The period is from 1min 19 hrs. 59 min

The time and speed are displayed separately on dual LCD screens

There are numerous inserts and fittings available

Designed with a brushless DC motor requires no maintenance

Specifications :

Ambient Temperature	5 to 40°C
Amplitude	20mm
Humidity	0.8
Maximum Load	3kg
Motor Type	DC motor
Speed Display	LCD
Speed Range	40 to 200rpm
Time setting range	0 to 19h 59min
Voltage	100 to 240V, 50/60Hz



Email: info@labdel.com | Website: www.labdel.com