



Water Steam Distiller TRWD-605

www.labdel.com | info@labdel.com

Water Steam Distiller TRWD-605

Water Steam Distiller TRWD-605 is a corrosion-resistant as well as aging-resistant stainless steel mesh tower mechanism. Distillation of the generating equipment, including the condenser, evaporator, electric heating pipe, and control unit. It has an intelligent, automatic control system. It comes with a control panel and a button-type mode. Designed with a water level window that displays the residual water level.

Features:

It is a water tower mechanism made from steel

It is corrosion-resistant and age resistant

It has an automatic and intelligent control system

designed with a water level window that displays the residual water level

The heating component is made up of immersion heating tubes that are long-lasting

Specifications :

Output	10 L/hour
Package	610×310×570mm
Power	AC380V/ 7.5 KW
Weight	11 Kg



Email: info@labdel.com | Website: www.labdel.com