



Fluorescence Spectrophotometer TRFS-601

www.labdel.com | info@labdel.com

Fluorescence Spectrophotometer TRFS-601

Fluorescence Spectrophotometer TRFS-601 has a luminous and fluorescence concentration that can be used for kinetic determination and quantity analysis. It has a wavelength of 365 nm and a precision of 1 nm. The C-T diffraction grating gives information such as excitation source spectroscopy within emissions monochromaticity of 200 to 900 nm.

Features:

There are two operation modes available fluorescence intensity and luminous intensity. Scanners for fluorescence Concentration, direct reading, auto 0 adjustment, auto background subtraction, and other features are available

? The graphic calculation, as well as storage access

? In fluorescence intensity mode, kinetic determination and quantity analysis could be performed

? We could regulate data in online mode by using amount and quality software

There are offline and online modes available. in offline mode

The fluorescence emission measurement is supplied by the device's computer system

As well as analysis via the USB2.0, interface

semi-automatic sample introduction accessories, membrane sample accessories, microscale capillary sample holder

Fluorescence sample holder with various features, 200l micro-scale centrifuge tube

Provide a PC qualitative software package with time and wavelength aspects involved as an optional extra

Measurement accessories, such as a single-hole cell holder, are available as options

Accessories for powder samples, jacket samples, and so on

This technique can deliver information counting excitation and emission spectrum

Specifications :

Display Value	0.00 to 600.0
Emission Monochromatic	C-T diffraction grating (EM 200 to 900 nm)
Exciting Optical Filters	Interface optical filter (10 nm bandwidth)
Interface	USB2.0
Linearity Deviation	±0.995
Power	220V±22v 50HZ±1HZ
Sensitivity	S/N≥150
Size	550×510×365 mm
Stability	Better than 1.5 %/10min
Wavelength	365 nm
Wavelength Accuracy	±1 nm
Wavelength Repeatability	±0.5 nm
Weight	13 kg



Email: info@labdel.com | Website: www.labdel.com